



Iron Duck
2012 / 2013 Catalog



THE WORLDS *BEST* EMERGENCY EQUIPMENT





From the President:

Iron Duck products are designed to meet very demanding standards - YOURS. Because in the field, there is no room for error, no award for second best. Only the highest quality products can survive the constant abuse of extreme environments - from the mean streets of our major cities to the harsh realities of wilderness search and rescue. That's why Iron Duck bags and kits are made from the most durable fabrics available, triple stitched for reinforcement, and finished with all metal hardware.

Task-specific design, superior materials and workmanship. Fast, error-free deployment in extreme environments. Iron Duck products are Built for Life and backed by our Iron Clad Warranty.

- Michael J. Fleming President & CEO

Mission Statement:

Iron Duck will be the world's leading provider of products to the emergency medical, fire fighting and law enforcement professional. This will be achieved through our commitment to fulfilling the needs of our customers, producing the finest products, and providing exceptional service in an environment and culture conducive to the enrichment of our company, employees and community for today and tomorrow.

To Our Valued Customers:

Iron Duck's headquarters is adjacent to Westover Air Force Base in Chicopee, Massachusetts. This location provides input to the requirements of the individual rescuer and large response teams. At Iron Duck we know that you need equipment that will perform and remain durable in the harsh environments that life saving rescues can take place in. Feel confident and proud that you are buying an American product crafted of the finest materials and workmanship anywhere.

- Tim Rogers National Sales Manager

Contact Iron Duck:

Iron Duck
A Division of Fleming Industries
20 Veterans Drive
Chicopee, MA 01022
United States of America

Toll Free: 1.800.669.6900
Main Line: 413.593.3300
Fax: 413.593.5800
Toll Free Fax: 1.800.757.8585
www.ironduck.com/contact



Product Warranty:

At Iron Duck we take pride in the quality, durability, and performance of the emergency medical product that you count on everyday. Iron Duck products are warranted to be free from defects in material and workmanship when used for the purpose intended under normal conditions. Damage caused by accident, abuse, improper care, and wear and tear is not covered by the warranty.

Product Care:

Treat all 1000D Nylon, Liner and Webbing components with a nylon safe fabric protector prior to first use and after each washing. Treat all zippers with the same fabric protector to ensure long life and smooth operation. Clean Universal Precautions fabrics with all-purpose soap and water.

**Do not use chlorine bleach.

**DO NOT Machine Wash Dirt or soil can be removed by wiping or brushing with household detergents (including antibacterial detergents) or cleaners.

Grease and oil stains can be removed with a commercial dry cleaning fluid. To clean a carry case contaminated by blood; Fresh stains can usually be removed by a thorough cleaning with cold* water and detergent. (*Hot water will "set" the stain) Persistent stains: Soak in a solution of 1 qt. warm water, 1.5 tsp dish washing or liquid laundry detergent, and 1 tb ammonia for 15 minutes. Brush or scrape gently to lift excess solid or caked-on stains, blotting occasionally with an absorbent pad, until the stain is removed. Rinse well with water, making sure to remove all traces of the ammonia. Wipe all exterior and interior surfaces with household detergent (antibacterial detergents can be used).

Cleaning:

To clean and disinfect a backboard contaminated by blood or other bodily fluids:

Wash and scrub the stain with a cleaning solution recommended by your local protocol. (Anti-bacterial solutions, Hydrogen Peroxide and Bleach in a 1:10 ratio are all products approved for use with Iron Duck backboards.) If stain is still visible, cover with solution and lay dry heavy-duty paper towels or brown paper (such as that used in the brown paper grocery bags) over the stain and cleaning solution, press down lightly to ensure that the solution permeates the paper.

Allow the board to dry overnight. Discard the used paper towels/wipes in biohazard bags.

The stain may require additional scrubbing with a stiff brush to remove debris left over from the soaking procedure.

Product Customization : Screen Printing for Bag & Board Customization:



Custom screen prints of simple text graphics in black ink are available on any Immobilization Backboard at a cost of \$10.00 per unit. If the custom screen print includes a logo, the digital artwork requirements are:

1. Black and white digital artwork (no gray scale; no lines under 1 pt thick) in a usable file format. These include:

a. Vector based file types are ideal since they can be resized with no loss of quality.

Vector based file types include .eps and .ai files.(Note: All text must be converted to outlines or the image flattened prior to sending.)

b. Pixel based file types are also acceptable as long as they have been saved at 300dpi at the approximate print size of 3" to 4". Pixel based file types include .pdf, .jpeg or .jpg, .tiff, and .psd files.

c. Logos saved from a website are usually .gif files that are 72 dpi and these file types cannot be used for screen printing.

2. Logos not meeting these requirements can be recreated for an additional fee of \$75.00 per hour. Fax or email the available file type to receive a quote.

Bag Embroidery Customization:



Custom screen prints of simple text graphics in white ink are available on any Emergency Medical Carry Cases at a cost of \$10.00 per unit. If the custom screen print includes a logo, the digital artwork requirements are:

1. Black and white digital artwork (no gray scale; no lines under 1 pt thick) in a usable file format.
 - a. Vector based file types are ideal since they can be resized with no loss of quality. Vector based file types include .eps and .ai files.(Note: All text must be converted to outlines or the image flattened prior to sending.)
 - b. Pixel based file types are also acceptable as long as they have been saved at 300dpi at the approximate print size of 3" to 4". Pixel based file types include .pdf, .jpeg or .jpg, .tiff, and .psd files.
 - c. Logos saved from a website are usually .gif files that are 72 dpi and these file types cannot be used for screen printing.
2. Logos not meeting these requirements can be recreated for an additional fee of \$75.00 per hour. Fax or email the available file type to receive a quote.

Spine Boards

Base Board - 35850



- > **Weight Capacity: 500 lbs.**
- > **Float Capacity: 250 lbs.**
- > **Dimensions: 72" x 16" x 1.75"**
- > **Weight: 16 lbs**
- > **100% X-Ray, MRI and CT Scan compatible**
- > **Rotomolded seamless HDPE shell over non-toxic Eco-Friendly foam**
- > **Fourteen over sized raised hand holds**

STANDARD COLORS:

Ultra Vue 16 - 35775



- > **Weight Cap: 1000 lbs.**
- > **Float Capacity 250 lbs.**
- > **Dimensions: 72" x 16" x 1.75"**
- > **Weight: 16 lbs.**
- > **100% X-Ray MRI and CT Scan compatible**
- > **Rotomolded one piece seamless HDPE shell**

STANDARD COLORS:

Ultra Vue 18 - 35755



- > **Weight Cap: 1000 lbs.**
- > **Float Cap: 250lbs.**
- > **Dimensions: 72" x 18" x 1.75"**
- > **Weight: 22 lbs.**
- > **100% X-Ray MRI and CT Scan compatible**
- > **Fourteen over sized raised hand holds**
- > **Available with or without 12 or 22 speed clip pins**
- > **Rotomolded 1 piece seamless HDPE shell over non-toxic Eco-Friendly foam**

STANDARD COLORS:

Pedi Air Align Complete - 35840



- > **Capacity: Pediatric Patients Aged 0-5 Years**
- > **Dimensions: 48" x 12" x 1.75"**
- > **Weight: 7.5lbs.**
- > **100% X-Ray, MRI and CT Scan compatible.**
- > **Designed with a patented Dual Plane Head Drop System.**

STANDARD COLORS:

Spine Boards

Ultra Short Board - 35940



- > Dimensions: 31.25" x 15.75" x 1"
- > Weight: 4 lbs
- > Design: Rotomolded one piece seamless HDPE shell over nontoxic Eco-Friendly foam
- > Six Over Sized Hand Holds
- > 4 Pediatric Strapping slots

STANDARD COLORS: Y B OR R LG OD BK

Ultra Spac- Sav - 35940



- > Weight Cap: 500 lbs. Weight Capacity
- > Dimensions: Standing: 72" x 16" x 2"
- > Folded: 36" x 16" x 4"
- > Weight: 19 lbs
- > Stores in half the space
- > Ideal for the long search and rescue, urban stair climb, patient transport vehicles and quads

STANDARD COLORS: Y B OR R LG OD BK

Ultra Loc - 35900



- > Weight Capacity: 1000 lbs.
- > Float Capacity: 250 lbs.
- > Dimensions: 72" x 16" x 1.75"
- > Weight: 14 lbs
- > 100% X-Ray, MRI and CT Scan compatible
- > Rotomolded seamless HDPE shell over non-toxic Eco-Friendly foam
- > Fourteen over sized raised hand holds
- > Available with or without 10 speed clip pins

STANDARD COLORS: Y B OR R LG OD BK

Immobilization

Board Loc | 30000



- > Fully Adjustable
- > Hook and Loop Closures
- > Includes a Nylon Pouch (Shown blue - color may vary)

Stretcher Straps



- Iron Duck straps are the most durable in the market. We offer this item in a variety of size, material and fasteners.
- Available as 5', 7', and 9, lengths
- 1 or 2 piece
- Plastic / Metal Buckle
- Loop or Speed Clip
- Nylon or UP

NYLON & UP COLORS: OR BK

Immobilization

Board Loc | 30000



Includes:

- > Carry Case: 14"L x 6"W x 8"H
- > 1 Black/Orange Pelvis Strap (78"L x 3"W)
- > 1 Black/Yellow Femur (73"L x 3"W)
- > 1 Black/Red Chest Strap (68"L x 3"W)
- > 1 Black/Green Upper Tibia Strap (63"L x 3"W)
- > 1 Black/Blue Lower Tibia Strap (58"L x 3"W)

STANDARD COLORS:  R

Board Loc | 30000



- > Fully Adjustable
- > Hook and Loop Closures
- > Includes a Nylon Pouch (Shown blue - color may vary)

STANDARD COLORS:  B

Disposable Adult Foam Head Blocks | 35979



- > The Iron Duck Disposable Head Block system features an adhesive back.
- Sold as 10 pack (2 blocks per pack)

STANDARD COLORS:

RB



Disposable Pedi Foam Head Blocks | 35845



- > The Iron Duck Disposable Head Block system features an adhesive back.
- Sold as 10 pack (2 blocks per pack)s

STANDARD COLORS:  RB

Ultra Space Save Case | 35940-CASE



- > Made of Universal Precautions material
- > Dimensions: 39"l x 16"w x 7"h
- > Standard color: Red (color may vary)

STANDARD COLORS:  R

IDEA | 35960



- > Standard extrication device procedures apply
- > Allows for spine and thorax immobilization
- > May also be used as a hip and pelvic splint
- > Replaceable head, body and leg straps
- > Body straps are color coded for instant identification

STANDARD COLORS:  R

EMS Backpacks

Ultra Backpack | 32440

STANDARD NYLON COLORS: RD  GN 



This is the largest backpack Iron Duck offers. It has exceptional organizational characteristics and is a durable multi-purpose pack. This workhorse is designed for quick deployment of AED, accessories, oxygen tank, adjuncts and trauma ALS supplies. Overall Dimensions: 23"H x 19"W x 12"D.

Ultra Backpack | 32440

STANDARD NYLON COLORS: RD  RB 



This is the midsize workhorse backpack of Iron Duck's line. This multi-purpose pack is designed for easy organization, secure transport, quick deployment, and houses an AED, accessories, oxygen tank, adjuncts, and trauma and ALS supplies. This unit boasts a slim profile, ergonomically designed carry options, and is made from our heavy duty nylon fabric. Overall Dimensions: 23"L x 14"W x 12"D.

BLS Backpack | 35132

STANDARD NYLON COLORS: RD  GN 



Our most compact Iron duck backpack offered. This multi-purpose case is designed very simply for easy organization, secure transport and quick deployment of an AED and accessories, oxygen tank and adjuncts, trauma and BLS supplies. A durable pack boasting a slim profile, multiple ergonomically designed carry options and an economical price. Overall Dimensions: 21"L x 14"W x 8"D.

OMNI II Master Case Backpack | 34025

STANDARD NYLON COLORS: OR  GN 







Iron Duck's Omni II Backpack is our most economical multi-purpose pack. It provides a simple yet complete organization of airway, O2 and trauma supplies in a versatile design. Overall Dimensions: 24"L x 13"W x 7.5"H.



EMS Trauma Bags




Ultra Sofbox Plus | 32325

STANDARD NYLON COLORS: RD  RB  OR 
STANDARD UP COLORS: RD 



This new style of the Sofbox Plus has tapered exterior pockets and a UP covered bottom with impact resistant runners for the ease of cleaning and provide a professional contemporary look. Featuring full length front and rear pockets that each house multiple interior pouches and elastic holders for supplies.
Dimensions: 23"L x 17"W x 12"H ~ 4692 cu in

Pack Case Triple | 32499AT

STANDARD NYLON COLORS: RD  NB  OR 
STANDARD UP COLORS: RD 



- > Dimensions: 19"L x 14.5"W x 8.5"H 2341.75 cu in
- > Rugged 1000D construction features polypro handles and a UP bottom.
- > A favorite "jump kit" for both BLS and ALS crews, ideal for first aid applications.



Pack Case Plus | 32499A

STANDARD NYLON COLORS: NB 



- › Includes shoulder strap
 - › Dimensions: 19"L x 14"W x 12"H ~ 3132 cu in
- his new style of the Trauma Pack Plus has tapered exterior pockets and a UP covered bottom for the ease of leaning and provide a professional contemporary look.

Trauma Pack Plus | 32350

STANDARD NYLON COLORS: RD  RB  OR 
STANDARD UP COLORS: RD 



- > Includes shoulder strap
 - > Dimensions: 19"L x 14"W x 12"H ~ 3132 cu in
- This new style of the Trauma Pack Plus has tapered exterior pockets and a UP covered bottom for the ease of cleaning and provide a professional contemporary look.

Trauma Bag Standard | 36001S

STANDARD NYLON COLORS: NB  OR 



- > Dimensions: 17"L x 9"W x 7"H ~ 1,040 cu in
- > Iron Duck's Standard has long been the industry standard for first responders. The full-length front pocket is split into two separate pouches and a row of elastic tacked at 2" intervals is positioned just above for easy access to small diagnostic and trauma tools.



EMS Oxygen Bags

Breathsaver | 34016D



- > Includes shoulder strap
- > Dimensions: 27"L x 13"W x 11"H ~ 3861 cu in

The Breathsaver is one of the most popular and versatile airway management bags in the market. This bag efficiently organizes tank, valve masks, airway adjuncts, hand-held suction and multiple trauma supplies.

STANDARD NYLON COLORS: RD  RB  OR  GN 
STANDARD UP COLORS: RD 

Breathsaver Plus | 34016DP



- > Includes shoulder strap
- > Access to oxygen tank from side pouch
- > Dimensions: 27"L x 13"W x 11"H ~ 3861 cu in

The Breathsaver Plus includes all of the features of the Breathsaver and more. This bag efficiently organizes tank, valve masks, airway adjuncts, hand-held suction and multiple trauma supplies.


STANDARD NYLON COLORS: RD  GN 
STANDARD UP COLORS: RD 

Ultra Breathsaver | 34018



- > Deluxe ergonomic dual adjustable shoulder strap included
- > Access to oxygen tank from side pouch
- > Dimensions: 29"L x 13"W x 12"H ~ 4524 cu in

The Ultra Breathsaver is the best trauma deployment system in the market. This bag was designed by field professionals for efficient storage and fast, error-free deployment of oxygen/resuscitation, airway management, and trauma supplies.

STANDARD NYLON COLORS: RD  RB  OR 
STANDARD UP COLORS: RD 

Oxygen Bag | 36002D



- > Includes shoulder strap
- > "D" Size 21"L x 8.5"W x 8.5"H ~ 1,517 cu in. This classically styled cylindrical clamshell design is ideal for protection and fast deployment of a D or Jumbo D tank, regulator and accessories. Fully lined with high impact foam the removable liner makes cleaning easy.

STANDARD NYLON COLORS: GN 
STANDARD UP COLORS: RD 

EMS Oxygen Bags

Oxygen Bag E-Size | 36002E

STANDARD NYLON COLORS: GN 
STANDARD UP COLORS: RD 



- > Includes shoulder strap
- > "E" Size 32"L x 8.5"W x 8.5"H ~ 2,312 cu in. This classically styled cylindrical clamshell design is ideal for protection and fast deployment of a E tank, regulator and accessories. Fully lined with high impact foam the removable liner makes cleaning easy.

Oxygen Bag 1-Pocket | 36002-PK

STANDARD NYLON COLORS: GN 
STANDARD UP COLORS: RD 



- > Includes shoulder strap
- > Dimensions: 21"L x 11.5"W x 8.5"H 1,989 cu in. This clamshell designed oxygen bag with 1 full length exterior pocket provides added storage and is ideal for protection and fast deployment of a D or Jumbo D tank, regulator and accessories.

Oxygen Bag 2-Pockets | 36002-2PK

STANDARD NYLON COLORS: GN 
STANDARD UP COLORS: RD 



- > Includes shoulder strap
 - > Dimensions: 21"L x 14"W x 8.5"H ~ 2,499 cu in
- This clamshell designed Oxygen bag with 2 full length exterior pocket provides added storage and is ideal for protection and fast deployment of a D or Jumbo D tank, regulator and accessories.

EMS Modules

Laryngoscope Bag | 33045

STANDARD NYLON COLORS: RB 



- > Dimensions: 14"L x 1.5"W x 10"H
- Another industry classic, Iron Duck's attache style zippered design holds Adult and Pedi Laryngoscope handles and a full range of blades. The attaché style zippered design consists of a durable 1000 denier nylon fully padded exterior with clear dividers inside that contain elastic tacking beneath to hold handles, blades, tubes and forceps securely without sacrificing easy identification.

Meds Pack | 32375

STANDARD NYLON COLORS: RB 



- > Dimensions: 13"L x 3"W x 5.5"H
- Iron Duck's Ultra Meds Pack is a compact drug module allowing separate storage, easy restocking and quick access to delicate ampoules and vials, whether boxed or unboxed. The exterior shell is constructed of 1000D Nylon and features a web quick-grab handle.

EMS Modules

IDX Expandable Drug Module | 34510

STANDARD NYLON COLORS: RD 



- > Standard Nylon Colors: Red
 - > Dimensions: 9"L x 6.5"W x 3.5"D (6"D fully expanded)
- Iron Duck's IDX Expandable Drug Module supports easy restocking and fast deployment via notebook style binding. The Drug Module holds up to 8 removable & customizable pages with multiple pockets and elastic tacking for supplies.

IDX Intubation Module | 34520

STANDARD NYLON COLORS: GN 



- > Dimensions: 13"L x 8"W x 2.5"H
 - > Dimensions open: 13"L x 27.5"W x 2.5"H
- Iron Duck's tri-fold Airway Module organizes a complete set of intubation supplies including laryngoscope handles and blades, airways, ET tubes and holders and accessories.

IDX IV Module | 34530

STANDARD NYLON COLORS: YL 



- > Dimensions: 7"L x 11"H x 1.5"D
 - > Open Dimensions: 21.5"L x 11"H x 1.5"D
- Iron Duck's tri-fold I.V. Module offers compact organization of a complete set of emergency I.V. supplies. The exterior shell is made of 1000Denier Nylon and features a quick-release plastic buckle closure and a clear vinyl card-sized window.

IV Starter Pack | 36015

STANDARD NYLON COLORS: RB 



- > Dimensions: 11"L x 7"W x 1.5"H
- The I.V. Starter Pack was designed by field professionals to hold all the necessary equipment for quick I.V. set-up. There is room for everything from catheters and tape, to blood drawing tubes and 10 cc syringes. Made from 1000 denier nylon and completely padded with closed cell foam and conveniently fits into all of our Oxygen and Trauma Cases.

IV Starter Caddy | 36017

STANDARD NYLON COLORS: RD 



- > Dimensions: 11"l x 4"w x 3.5"h
- Made from 1000 denier nylon this self contained unit designed for multi-casualty incidents. Holds up to a 1000ml I.V. bag with room to spare. Inside full-length pocket along side wall allows easy storage of I.V. tubing and blood drawing equipment. Six exterior pockets provide storage for catheters. End pocket is large enough to hold 1" tape and numerous alcohol preps, band-aids, etc.



EMS Midwife

Breathsaver-34016D-PR-Midwife

STANDARD NYLON COLORS: PR 



- > Includes shoulder strap
- > Dimensions: 27"L x 13"W x 11"H ~ 3861 cu in

The Breathsaver is one of the most popular and versatile airway management bags in the market. This bag efficiently organizes tank, valve masks, airway adjuncts, hand-held suction and multiple trauma supplies. The large main compartment contains tank straps and mesh pouches for organization and has a removable liner for easy cleaning. Heavy duty handles wrap around the entire bag for better weight distribution.

Trauma Pack Plus-32350-PR-Midwife

STANDARD NYLON COLORS: PR 



- > Includes shoulder strap
- > Dimensions: 19"L x 14"W x 12"H ~ 3132 cu in

This new style of the Trauma Pack Plus has tapered exterior pockets and a UP covered bottom for the ease of cleaning and provide a professional contemporary look. Featuring full length front and rear pockets, the front pocket has multiple interior pouches and elastic holders for supplies while the rear pocket has elastic holders and a generous volume of storage space.

Ultra Sofbox Plus-32325-PR-Midwife

STANDARD NYLON COLORS: PR 



- > Dimensions: 23"L x 17"W x 12"H ~ 4692 cu in

This new style of the Sofbox Plus has tapered exterior pockets and a UP covered bottom with impact resistant runners for the ease of cleaning and provide a professional contemporary look. Featuring full length front and rear pockets that each house multiple interior pouches and elastic holders for supplies.



www.ironduck.com / 1.800.669.6900